

Specimen Coating for High Resolution Electron Microscopy

Phil Missing, Emitech, Ltd.

Probably the most common coating technique in electron microscopy is sputter coating, primarily for SEM where a thin electrically conducting film is required to prevent charging and to enhance secondary electron emission.

These systems are relatively inexpensive, and the modern instruments are easy to use, often being fully automated.

D.C. Sputtering

The direct current (dc) sputter coater uses the arrangement of two electrodes. One positive - the specimen stage, and one negative - the target, typically gold. A glow discharge is formed, and the sputtered material from the target is deposited onto the specimen.


The sputtering efficiency is often enhanced by the use of a magnetron head. An inert gas such as Argon is used, and a suitable low pressure is achieved using a two stage rotary vacuum pump. The grain size for gold target material is of the order of 2 - 5 nm. To achieve a continuous film, the thickness of deposition is typically 10 nm but may be somewhat thicker, up to 20 nm for an irregular specimen.

Chromium Deposition

With a turbo sputtering system it is now possible to sputter chromium and achieve fine grain size of order of 0.5 nm thick - which is electrically conductive and continuous. This development in sputter coaters provides ultra thin, high quality films for the users of high resolution SEM.

Until recently, these coatings have been suitable for most SEM applications. With the advent of high resolution SEM, in particular field emission, there is an increasing requirement to improve the performance of the coatings to smaller, continuous films.

While ion beam systems are being employed, the deposition times are long and the instrumentation is relatively expensive.




Ultra Thin Continuous Chromium Film

Turbo Sputtering

There are benefits in depositing ultra thin films of chromium on SEM specimens to give improved resolutions at higher magnifications.

To avoid the formation of oxide and contamination when using chromium as the target material, it is necessary to operate at pressures lower than those achievable using two stage rotary vacuum pumps. Such systems usually employ a turbomolecular pumping system.

The system is pumped down to 1×10^{-5} mbar. Using an argon bleed, sputtering is typically at 1×10^{-2} mbar. ■



Emitech K575 Chamber Arrangement

ANALOG/DIGITAL PRINTERS



Seikoshia Video Printers

- 256 Grey Levels
- 10 MHz - 135 MHz
- 5¢ - 19¢ - 40¢ Hardcopies
- 3 x 4, 8 1/2 x 8 or 8 1/2 x 11 Formats



Mitsubishi Printers

- Dye-Sublimation Printers
- Analog/Digital Input
- 16.7 Million Colors
- 4 x 4 thru 12 x 18 Format
- Thermal Greyscale Printers



Tektronix Printers


- Dye-Sublimation Printers
- Digital/Postscript Input
- 16.7 Million Colors
- 8 1/2 x 11 Thru 8 1/2 x 14
- TCP/IP and Ethertalk

VITAL IMAGE TECHNOLOGY INC. OFFERS COMPLETE SALES AND SERVICE, INCLUDING FREE INSTALLATION IN MOST PARTS OF THE UNITED STATES.



VIT

(800) 860-IMAGE



VITAL IMAGE TECHNOLOGY INC.

MICRO STAR DIAMOND KNIVES

★ BEST QUALITY

Backed the best guarantee: You don't pay until you have tested your knife and are totally satisfied with its flawless quality.

★ BEST SERVICES

Exchange (or sharpen) your old knife, any brand, type, size or age for a new MICRO STAR at the sharpening price.

★ BEST SELECTION

8 boat styles, 15 types and 59 sizes, from 1 to 10mm. The most popular sizes: 2, 3, 4 and 5 mm ready for



IMMEDIATE DELIVERY


800 533 2509 FAX: 409 294 9861 E-MAIL: US3SNQ7N@IBMMAIL.COM

EMITECH K SERIES


THE COMPLETE 'SPECTRUM' OF PREPARATION EQUIPMENT



*K575 (illustrated) Turbo Pumped Sputter Coater for Chromium.
K550 Sputter Coater for Gold.
K650 Large Sample Sputter Coater.*




*K775 (illustrated) Turbo Pumped Low Temperature Freeze Drier.
K750 Peltier Stage Freeze Drier.*



K850 Critical Point Drier.



*K975 (illustrated) Turbo Pumped Evaporator System.
K950 Turbo Pumped Carbon Evaporator.
K450 Carbon Fibre Coater.*



K1250 Cryogenic System with Integral Turbo Pumping, and Interfacing to S.E.M.

SERVING THE SCIENCE OF ELECTRON MICROSCOPY

For further details please contact:
LINDA DAILEY, EMITECH USA, 3845 FM 1960-WEST, SUITE 345, HOUSTON, TEXAS 77068 1-800-444-3137 (713) 893 2067 FAX: (713) 893 8443
BOB SNELL, CANBERRA-PACKARD CANADA Ltd. 6470 VAN DEEMTER COURT, MISSISSAUGA, ONTARIO L6T 1S1 1-800-387-9559 (416) 795 2559 FAX: (416) 795 2598

Circle Reader Inquiry #10