



CORPORATE PARTNER

WHITE PAPER

## Non-destructive Hall measurement environment facilitates novel device R&D

LAKE SHORE CRYOTRONICS, INC.

Andy Phillips

Hall measurement systems are invaluable research tools for electronic materials and device characterization, providing quantitative data on carrier type, concentration, and mobility, as well as measurements related to resistivity,  $I$ - $V$  curves, and other useful information. Understanding the underlying transport mechanisms of a material provides insights into its potential for use in new device applications.

Sample preparation for Hall measurements can be a time-consuming process. Conventional Hall measurement techniques generally require a small, dedicated sample with hard-wired electrical connections. This typically entails destruction of at least part of the sample and prevents repeated measurement of the same structure at different stages of device fabrication. In addition, while Hall systems generally offer measurement at a range of temperatures, their sample environments may not be ideally suited to materials that react quickly in atmosphere. In these cases, a different measurement platform is indicated.

Recognizing the need for non-destructive testing of wafer-scale materials, and with an eye toward providing improved environmental control for Hall measurements, Lake Shore has been working to integrate its vertical field-equipped cryogenic probe station platform with its Hall measurement instrumentation. This type of configuration makes for easier sample preparation and enables measuring of multiple structures in a single experimental session, as well as providing the ability


to repeat measurements of Hall structures at progressive steps in the formation of devices co-located on the wafer.

Since carrier densities and transport mechanisms change with temperature in very specific ways, a range of Hall measurements starting at cryogenic temperatures are typically taken. Cryogen-free refrigeration and advancements in thermally compliant cryogenic probe design now enable automated variable temperature measurement runs while ensuring consistent contact quality. The probe station's tightly controlled vacuum environment accommodates materials that are reactive to air and those that require initial warming to drive out moisture. Finally, the probe station platform allows the researcher to easily reconfigure the attached instrumentation in order to make any other (non-Hall) electrical measurements required while the sample is still in the environment.

### Easy sample handling

In conventional Hall measurement systems, sample preparation usually involves dicing a wafer into appropriately sized pieces and attaching wires by soldering, wire bonding, or, if all else fails, conductive epoxies. Soldering and wire bonding can be especially tricky with high-bandgap materials like SiC, while sensitive materials like organics can be damaged in the process. The use of a probe station platform eliminates the need to attach wires to the sample, significantly speeding sample prep and avoiding any damage to the sample. Instead, four or six micromanipulated probes are individually positioned as needed once the sample has been installed in the sample chamber.

The probe station provides a large sample stage that can accommodate multiple individual samples and is ideal for handling wafer-scale samples (in the case of the Lake Shore system, up to 51 mm in diameter). In the typical case where the device developer has patterned a variety of different test structures on the same wafer, the repositionable probes can shift from measuring one Hall structure to the next without affecting the overall sample environment. Transport characteristics can be quickly confirmed at the current stage of device development, and the wafer remains intact for subsequent processing steps (see Figure 1).



**Figure 1.** A probe station sample stage can hold multiple samples, and repositionable probes simplify the process of measuring multiple Hall structures.

Andy Phillips, Lake Shore Cryotronics, Inc. The MRS Corporate Partner Program supports the Materials Research Society Foundation.

### A controlled environment

The probe station environment has tightly controlled parameters of vacuum, temperature, and magnetic field. A vacuum chamber, capable of sustaining pressures of  $10^{-5}$  Torr and below, offers an advantage over conventional Hall measurement systems. Materials can be warmed in vacuum prior to measurement to drive out any residual moisture. Reactive materials can be measured in an environment that minimizes atmospheric exposure and resulting contaminations. Cryo-pumping techniques resulting from creating temperature differentials within the probe station drive contaminants from the sample space.

Variable temperature Hall measurement experiments are commonly used to reveal specific changes in transport characteristics. At cryogenic levels, carriers can be identified by their activation energies and provide clues to predominant transport mechanisms in materials. Hall measurement at these low temperatures is especially helpful in understanding the frequency response of the material or device. A cryogen-free refrigerator and cryogenic temperature controller ensure reliable closed-loop management of sample temperature over a broad range. For best thermal accuracy, probes and probe arms need to be thermally anchored to the sample stage, while use of radiation shielding reduces the amount of stray thermal energy impacting the sample.

Specialized, thermally compliant probes (such as Lake Shore's ZN50R-CVT) enable automated measurement over a wide temperature range. These probes are designed to maintain consistent contact with the material despite thermal expansion of the sample, eliminating the need to manually adjust position as experiment temperature rises and allowing for continuous variable temperature Hall measurements.

The high field superconducting magnet of the probe station (capable of around 2 T) facilitates advanced quantitative carrier analysis. For material with a mobility of  $5000 \text{ cm}^2/(\text{Vs})$  or higher, the presence of multiple carrier species in a sample can be determined. This is particularly useful


for resolving individual carrier mobilities and densities in quantum wells and high-electron-mobility transistors (HEMTs).

By employing specialized software, such as the Lake Shore QMSA package, properties of majority carriers and other carrier mechanisms of a multicarrier sample material can be accurately identified and analyzed. The software can segregate the mobility spectrum for each carrier species in multilayer structures, and is particularly useful where it is necessary to obtain magnetic field- and temperature-dependent information for the separate carriers within a material. Hall coefficient and resistivity can be measured at different magnetic field levels

and then correlated to a unique mobility spectrum. This multicarrier measurement analysis can be very helpful for troubleshooting designs, for instance, materials fabricated in R&D environments and to evaluate deposition products in quality control applications. Without a QMSA type of analysis, researchers are left with performing fixed-field Hall measurements that extract carrier properties in bulk, limiting their utility as an analysis tool.

### Example research application

Favorable results have been achieved in testing of HEMT devices. Considered highly promising for high-speed communications devices and applications, HEMT devices use a lattice-matched heterojunction between donor and acceptor compound semiconductors. Early HEMTs were largely fabricated using GaAs (gallium arsenide) substrates; advancing research is pursuing other semiconductor materials. In operation, the HEMT forms a two-dimensional electron gas (2DEG) layer that contains the mobile charge carriers and provides conduction for the device. With the aid of QMSA analysis, researchers have found a number of interesting properties




for these devices, including significant differences in carrier mobility with temperature change and unexpected differences in majority and minority carrier mobility. Hall coefficient and resistivity have been measured at different magnetic fields up to 2 T and then correlated to a unique mobility spectrum.


For its test application, Lake Shore used its Model 8425 DC Hall measurement system with integrated cryogenic probe station to measure 2DEG materials on a 2-inch (51-mm) AlGaIn/GaN patterned wafer (Figure 2). In this application, temperature-dependent behaviors were measured on three various sized device structures patterned on the same wafer, with Hall resistance as a function of applied field measured during one experimental session (no QMSA analysis was used for this 2DEG characterization). Figure 3 shows the results of the Hall resistance measurements over a range of fields for the devices. Blue data points show measurements at 10 K, while red data points show measurements at 305 K. The linear behavior of the Hall resistance versus field indicates that this material may have possible future use in sensing applications. □

## There's No Other AFM Like Cypher™

The Highest Resolution Fast Scanning AFM



Atomic point defects on calcite crystal face



Growth at a screw dislocation in a calcite crystal

Elastic modulus of a polymer blend

- Routinely achieves higher resolution than other AFMs
- The only full-featured fast scanning AFM
- blueDrive™ photothermal excitation reinvents tapping mode
- Effective, hassle-free environmental control

See us in booth #510

[www.AsylumResearch.com/NoOtherAFMLikeIt](http://www.AsylumResearch.com/NoOtherAFMLikeIt)



*The Business of Science®*

+1 (805) 696-6466  
[sales@AsylumResearch.com](mailto:sales@AsylumResearch.com)  
[www.AsylumResearch.com](http://www.AsylumResearch.com)



give to the  
causes you  
care about



## Education & Outreach. University Chapters. Awards.

### Your donations to the **Materials Research Society Foundation** have...

- contributed to the SciBridge where experiment kits connect African university students to US materials scientists
- supported "Science Saturdays," a series of science workshops for high school students in the Philadelphia area
- established a new MRS Award recognizing postdoctoral scholars
- provided education modules on sustainable materials for elementary and middle school students with diverse ethnic, social and economic populations
- brought hands-on science lessons to rural middle schools in the state of Tennessee through the Materials Outreach for Rural Education (MORE) program
- provided free MRS memberships for university students in developing countries
- honored outstanding contributors to the progress of materials research through 12 annual MRS awards
- expanded the MRS University Chapters Program globally, now with 89 Chapters around the world

# Become an MRS Congressional Science and Engineering Fellow!

Help **improve the interface** between science and legislative decision-making.

**Advocate for policies** that will facilitate the discoveries of the future.

Play a crucial role as you **educate the public** about the benefits of science.

Decisions made by Congress, regulatory agencies and local government have profound effects on the way in which science is conducted. By keeping decision-makers well informed on the current affairs of the scientific community, MRS Congressional Science and Engineering Fellows ensure the right choices are being made. Now's your time to make a difference!

The Materials Research Society offers material scientists two exciting opportunities to participate in, and contribute to, the federal policymaking process, while learning firsthand about the intersection of science and policy.

During your year as a Fellow you will:

- ▶ contribute widely to the effective use of materials science knowledge in government
- ▶ broaden awareness about the value of scientist- and engineer-government interaction among society members and within government
- ▶ have significant freedom to follow specific topics and issues that interest you

The MRS Congressional Science and Engineering Fellowship Program is an invaluable experience, but don't just take our word for it. Our past Congressional Fellows explain it best!

"Academia taught me how to think, but the MRS Congressional Fellowship taught me how to get things done. Never have I had such leverage, such opportunities to comingle with dignitaries, to structure agreements and broker deals, as I did in that year. I learned how to navigate past armies of secretaries shielding a VIP, enlist military support for a project, take a rough idea and make it law, to fashion an event into a sound bite and then watch it propagate across the news. I learned to take data and present it in such a way that it gravitated, almost of its own accord, all the way up to the Vice President of the United States. These are skills anyone, who is going anywhere, can use."

**Merrilea Mayo** Founder, Mayo Enterprises, LLC.  
MRS Congressional Fellow 1998-1999  
Office of Senator Lieberman

"At the end of the fellowship year I found that I was enjoying "doing" science policy more than just teaching about it, and I ended up staying on in Rep. Honda's office as a member of the staff where I have remained for over a decade. I would not have had that opportunity without the Congressional Fellowship. I encourage anyone who wonders about how federal policies are developed or wants to have a greater role in that process to apply to be a Congressional Fellow."

**Eric Werwa** Legislative Director, Congressman Mike Honda  
MRS Congressional Fellow 2001-2002  
Office of Congressman Mike Honda

**MRS**® To learn more about the MRS Congressional Science and Engineering Fellowship Program and how you can apply, visit [www.mrs.org/congressional-fellows](http://www.mrs.org/congressional-fellows).

Applications for the 2015-2016 MRS Congressional Science and Engineering Fellowship Program are posted on the MRS website. **Deadline for submission is January 2, 2015.**