



Powder Diffraction

An international journal of materials characterization

Volume 15 Number 1 March 2000

33-1161
SiO₂
Silicon Oxide
Quartz, syn


Rad. CuKα₁ λ 1.40598 Filter Mono. d-sp Diff.
Cut off Int. Diffractometer $\theta/\lambda_{\text{Cu}}$ 3.6
Ref. Natl. Bur. Stand. (U.S.) Monogr. 25, 18 61 (1981)


Sys. Hexagonal
a 4.9133(2) b c 5.4053(4)
α β
Ref. Ibid. S.G. P3₂21 (154)
A Z 3 mp C 1.1001

D₁ 2.65 D₂ 2.66
Ref. Swanson, Fuyat, Natl. Bur. Stand. (U.S.) Monogr. 25, 18 61 (1981)

Color Colorless
Pattern taken at 25°C. Sample from the Glass Section at NBS.
Gathersburg, Maryland, USA. Ground single-crystals of optical qual-
ity. Pattern reviewed by Holzer, J., McCarthy, G., North Dakota State
University, Fargo, North Dakota, USA. ICDD Grant-in-Aid (1990).
Agrees well with experimental and calculated patterns. O₂Si type.
Quartz group. Also called: silica. Also called: low quartz. Silicon used
as internal standard. PSC: hP9. To replace 5-490 and validated by cal-
culated pattern. Plus 6 additional reflections to 0.9089.

dÅ	Int	hkl	dÅ	Int	hkl
4.257	22	100	1.1532	1	311
3.342	100	101	1.1405	<1	204
2.457	8	110	1.1143	<1	303
2.282	8	102	1.0813	2	312
2.237	4	111	1.0635	<1	400
2.127	6	200	1.0476	1	400
1.9792	4	201	1.0438	1	105
1.8179	14	112	1.0347	<1	401
1.8021	<1	003	1.0150	<1	214
1.6719	4	202	0.9898	1	223
1.6591	2	103	0.9873	1	402
1.6082	<1	210	0.9783	1	313
1.5418	9	211	0.9762	<1	304
1.4536	1	113	0.9636	1	320
1.4189	<1	300			
1.3820	6	212			
1.3752	7	203			
1.3718	8	301			
1.2880	2	104			
1.2558	2	302			
1.2285	1	220			
1.1999	2	213			
1.1978	1	221			
1.1843	3	114			
1.1804	3	313			





Microdiffraction

Texture

- Real-time data collection with unique 2-D HI-STAR detector
- Lowest background level for highest signal sensitivity
- Simultaneous large χ and 2θ data acquisition for fastest powder and texture investigations
- 2-D texture, residual stress, and microdiffraction with incomparable quality
- Small Angle X-ray Scattering capability for 2-D large structure examination
- Cross-coupled Göbel Mirrors for highly intense, point like parallel beam
- Comprehensive range of stages for every kind of sample
- Sophisticated laser/video system for easy sample alignment and positioning
- GADDS and DIFFRAC^{plus} software suite for every kind of application
- Unrivalled flexibility and versatility inherent in the D8 diffraction family

*Preferred Orientation,
Microdiffraction,
Grained Materials,
Large Structures...*

*...perfect for
D8 DISCOVER
with GADDS*

**find out
what's inside**



BRUKER ADVANCED X-RAY SOLUTIONS

Worldwide: Tel. (+49) 721/595-2888 • Fax (+49) 721/595-4587 • <http://www.bruker-axs.de> • E-mail: info@bruker-axs.de
 USA and Canada: Tel: (+1) 608/276-3000 • Fax: (+1) 608/276-3006 • <http://www.bruker-axs.com> • E-mail: info@bruker-axs.com

crystallographica searchmatch

The makers of *Crystallographica* are proud to announce the launch of *Crystallographica Search-Match*, an all-new search-match program for Windows 95/98/NT.

- ◆ Works with all versions of the Powder Diffraction File including the new cPDF
- ◆ Search using full powder diffraction pattern and/or peak list
- ◆ Automatic residual search for multi-phase identification
- ◆ Unique integrated Boolean card retrieval and display
- ◆ Single / multi-phase full pattern powder simulations
- ◆ Reads common file formats
- ◆ Built-in tools include peak finding and background / $k\text{-}\alpha_2$ stripping
- ◆ Report writing directly to Microsoft Word
- ◆ Full 32-bit technology delivering unrivalled speed and power

**Visit our new web site for details
or contact us for a free demo CD!**

Search-Match?

Upgrade to the state-of-the-art




OxfordCryosystems

3 Blenheim Office Park, Lower Road, Long Hanborough
Oxford · OX8 8LN · UK
Tel: +44 (0)1993 883488 · Fax: +44 (0)1993 883988
E.mail: info@OxfordCryosystems.co.uk

www.crystallographica.co.uk

X-ray diffraction for today and tomorrow



The inevitability of change has long been recognized - and the pace of change is accelerating at an ever increasing rate. This is true in the laboratory domain as in every other aspect of modern life. What you require from your X-ray diffractometer could be a world apart from your expectations today.

The Philips X'Pert PRO diffractometer system is an X-ray diffraction system with a virtually limitless capacity to adapt to changing circumstances. The PreFIX (Pre-aligned, Fast-Interchangeable X-Ray optics) system enables instruments to be reconfigured in a few minutes to handle different types of analysis.

A totally open design architecture guarantees full freedom to accommodate further advances in years to come.

For more information:
Philips Analytical
Lelyweg 1,7602 EA,
Almelo, The Netherlands
Tel. : +31 (546) 534444
Fax : +31 (546) 534592

www.analytical.philips.com

Denver X-ray Conference Booth #56, 57



PHILIPS

Let's make things better

Published online by Cambridge University Press

Editor-in-Chief

Ting C. Huang
6584 Radko Drive
San Jose, CA 95119-1924 U.S.A.
ting_huang@ibm.net

Managing Editor

Shannon Mataboni
JCPDS-International Centre for Diffraction Data
12 Campus Blvd.
Newtown Square, PA 19073-3273 U.S.A.
mataboni@icdd.com

Editor for Reviews and Reprints

Deane K. Smith
1652 Princeton Drive
State College, PA 16803-3273 U.S.A.
smithdeane@aol.com

Editor for New Diffraction Data

William E. Mayo
Rutgers University
Ceramics Department
Piscataway, NJ 08855-0909 U.S.A.
edisonjam@aol.com

Western European Editor

Norberto Masciocchi
Università di Milano
Dipartimento di Chimica
Strutturale e Stereochimica Inorganica
Via Venezian, 21
20133 Milano, Italy
norbert@csmtbo.mi.cnr.it

Eastern European Editor

Jaroslav Fiala
SKODA Research Ltd.
31600 Plzeň, Czech Republic
jaroslav.fiala@vsb.cz

Editor for Southeast Asia

Brian H. O'Connor
Curtin University
GPO Box U 1987, Perth 6001
Western Australia, Australia
toconnorb@cc.curtin.edu.au

Editor for Japan

Hideo Toraya
Ceramics Research Lab
Nagoya Institute of Technology
Asahigaoka, Tajimi 507 Japan
toraya@crl.nitech.ac.jp

International Reports Editor

Winnie Wong-Ng
National Institute of Standards and Technology
100 Bureau Drive Stop 8520
Gaithersburg, MD 20899-8520 U.S.A.
winnie.wong-ng@nist.gov

Editorial Advisory Board

P. Bayliss, Sydney, Australia
C. Z. Bojarski, Katowice, Poland
A. Brown, Dorset, England
D. Cox, Upton, New York
W. Eysel, Heidelberg, Germany
L. Frevel, Midland, Michigan
P. Gado, Budapest, Hungary
H. Goebel, Munchen, Germany
R. Jenkins, Newtown Square, Pennsylvania
G. G. Johnson Jr., State College, Pennsylvania
Q. Johnson, Livermore, California
J. I. Langford, Birmingham, U.K.
D. Louër, Rennes, France
G. J. McCarthy, Fargo, North Dakota
H. F. McMurdie, Gaithersburg, Maryland
M. E. Mrose, Gaithersburg, Maryland
M. Nichols, Livermore, California
R. L. Snyder, Columbus, Ohio
T. Yamanaka, Tokyo, Japan
R. A. Young, Atlanta, Georgia

AIP Production: Lin Miller, *Editorial Supervisor*
Andrea Witt, *Journal Coordinator*
Kelly Quigley, *Chief Production Editor*

Powder Diffraction is a quarterly journal published by the JCPDS-International Centre for Diffraction Data through the American Institute of Physics (AIP). *Powder Diffraction* is a journal of practical technique, publishing articles relating to the widest range of application—from mineral analysis to epitaxial growth of thin films and to the latest advances in software. Although practice will be emphasized, theory will not be neglected, especially as its discussion will relate to better understanding of technique.

Submit manuscripts (3 copies) to the most appropriate *Powder Diffraction* Editor listed on this page. The Editors will consider all manuscripts received, but assume no responsibility regarding them. Materials will be returned only when accompanied by appropriate postage. There is no publication charge. See *Powder Diffraction Notes for Authors* for additional information.

Proofs and all correspondence concerning papers in the process of publication should be addressed to: Editorial Supervisor, *Powder Diffraction*, AIP, Suite 1N01, 2 Huntington Quadrangle, Melville, NY 11747-4502.

For advertising rates and schedules contact AIP Advertising Department. Orders, advertising copy, and offset negatives should be sent to: Advertising Department, American Institute of Physics, Suite 1N01, 2 Huntington Quadrangle, Melville, NY 11747-4502; phone: (516) 576-2440; fax: (516) 576-2481.

Subscription Prices (2000)

	U.S.A. & Canada	Mexico, Central & South America	Europe, Mid-East & Africa*	Asia & Oceania*
Individual	\$60	\$85	\$85	\$85
Institutional or Library	\$105	\$105	\$105	\$105

*Subscription rates to Eastern Hemisphere include air freight service.

Back-Number Prices. 2000 single copies: \$30. Prior to 2000 single copies: \$30.

Subscription, renewals, and address changes should be addressed to *AIP Circulation and Fulfillment Division (CFD)*, Suite 1N01, 2 Huntington Quadrangle, Melville, NY 11747-4502. Allow at least six weeks advance notice. For address changes please send both old and new addresses and, if possible, include a mailing label from the wrapper of a recent issue.

Claims, Single Copy Replacement and Back Volumes: Missing issue requests will be honored only if received within six months of publication date (nine months for Australia and Asia). Single copies of a journal may be ordered and back volumes are available in print or microform. Individual subscribers please contact AIP Circulation and Fulfillment Division (CFD) at (516) 576-2288; (800) 344-6901. Institutional or library subscribers please contact AIP Subscriber Services at (516) 576-2270; (800) 344-6902.

Reprint Billing: Contact: AIP Circulation and Fulfillment Division, Melville, NY 11747-4502; (516) 576-2230; (800) 344-6909.

Copying: Single copies of individual articles may be made for private use or research. Authorization is given (as indicated by the Item Fee Code for this publication) to copy articles beyond the use permitted by Sections 107 and 108 of the U.S. Copyright Law, provided the copying fee of \$6 per copy per article is paid to the Copyright Clearance Center, 222 Rosewood Drive, Danvers, MA 01923, USA. Persons desiring to photocopy materials for classroom use should contact the CCC Academic Permissions Service. The Item Fee Code for this publication is 0885-7156/2000 \$6.00.

Authorization does not extend to systematic or multiple reproduction, to copying for promotional purposes, to electronic storage or distribution, or to republication in any form. In all such cases, specific written permission from AIP must be obtained.

Permission for Other Use: Permission is granted to quote from the journal with the customary acknowledgment of the source. To reprint a figure, table, or other excerpt requires the consent of one of the authors and notification to AIP.

Requests for Permission: Address requests to AIP Office of Rights and Permissions, Suite 1N01, 2 Huntington Quadrangle, Melville, NY 11747-4502; Fax: 516-576-2450; Telephone: 516-576-2268; E-mail: rights@aip.org.

Document Delivery: Copies of journal articles can be ordered for online delivery from the new Articles in Physics online document delivery service (URL: <http://www.aip.org/articles/>).

Reprints: Reprints can be ordered with or without covers only in multiples of 50 (with a minimum of 100 in each category) from AIP, Circulation & Fulfillment/Reprints, Suite 1N01, 2 Huntington Quadrangle, Melville, NY 11747-4502; Fax: 516-349-9704; Telephone: 800-344-6909 (in U.S. and Canada), or 516-576-2234.


Powder Diffraction (ISSN: 0885-7156) is published quarterly (4X annually) by the JCPDS-International Centre for Diffraction Data through the American Institute of Physics. 1999 Subscription rates: US\$105. POSTMASTER: Send address changes to *Powder Diffraction*, AIP Circulation & Fulfillment Division, Suite 1N01, 2 Huntington Quadrangle, Melville, NY 11747-4502. Periodicals postage paid at Huntington Station, NY 11746, and additional mailing offices.

Online Availability: Abstracts of journal articles published by the AIP and Member Societies (and several other physics publishers) are available from AIP's SPIN database, via AIP's Online Journal Publishing Service (OJPS) (<http://ojps.aip.org>).

Copyright © 2000 JCPDS-International Centre for Diffraction Data, 12 Campus Blvd., Newtown Square, PA 19073-3273. All rights reserved.
www.icdd.com/products/journals.htm

D/MAX-

2000/PC



- **Standard features provide the most flexibility for your investment**
- **Flexible optical system**
- **User-friendly software and hardware**
- **Attachments for every application**

Rigaku's D/MAX family of x-ray diffractometers continues to evolve with the 2000/PC.

The D/MAX-2000/PC includes focusing and parallel-beam optics, which can be exchanged by the user in half a minute without realignment.

Variable slits are standard, providing exceptional low-angle performance down to 1.5°, a real benefit for many pharmaceutical and geological applications.

Rigaku's graphite diffracted-beam monochromator is specially designed for both focusing and parallel-beam configurations — no need for additional hardware and time-consuming reconfiguration.

In addition, the D/MAX-2000/PC offers the highest degree of automation found in any diffraction system.

Whether you require sealed-tube sources, high-reliability rotating anodes, or mirror optics, Rigaku has it all.



North America:
9009 New Trails Drive
The Woodlands,
Texas 77381 USA
Tel: 281-367-2685
Fax: 281-364-3628
e-mail: info@rigaku.com
<http://www.rigaku.com/>

Japan:
3-9-12 Matsubara-cho, Akishima-shi,
Tokyo 196-8666, Japan
Tel: 81-42-545-8189
Fax: 81-42-545-7985
e-mail: rinttyo@rigaku.co.jp
Telex: J26418 RINT

Horst Joseph Pentinghaus, Michael Tovar	Obituary: Walter Hans Eysel·1935-1999	1
--	---------------------------------------	---

TECHNICAL ARTICLES

Rahul Surana, Raj Suryanarayanan	Quantitation of crystallinity in substantially amorphous pharmaceuticals and study of crystallization kinetics by X-ray powder diffractometry	2
J. C. Taylor, I. Hinczak, and C. E. Matulis	Rietveld full-profile quantification of Portland cement clinker: The importance of including a full crystallography of the major phase polymorphs	7
Hiroshi Hayakawa, Toshimi Takada, Hirotohi Enoki, and Etsuo Akiba	Crystallographic data of new superlattice phases for spinel LiMn_2O_4 at low temperatures	19
A.-E. Gobichon, J.-P. Auffrédic, and D. Louër	A temperature-dependent powder diffraction study of chromium lanthanum nitrate, $\text{LaCr}(\text{NO}_3)_6 \cdot 12\text{H}_2\text{O}$	23
L. De Stefano, F. De Luca, G. Buccolieri, and P. Plescia	Milling effects upon quantitative determinations of chrysotile asbestos by the reference intensity ratio method	26
Martin Oetzel, Franz-Dieter Scherberich, and Gernot Heger	Heating device for high temperature X-ray powder diffraction studies under controlled water vapour pressure (0–1000 mbar) and gas temperature (20–200 °C)	30
Wei Lianhu, Xiao Yunde, Zhang Jinbei, and Lin ShaoFan	An introduction of chemical structure coding in the X-ray Powder Diffraction File	38
Fred J. Wicks	Status of the reference X-ray powder-diffraction patterns for the serpentine minerals in the PDF database—1997	42


NEW DIFFRACTION DATA

A. Bentayeb, O. Masson, Y. Noack, A. Yacoubi, and A. Nadiri	A mineralogical study of a natural pyrophyllite	51
Umberto Costa, Paolo Ballirano	Improved powder X-ray data for the cement phase $\text{Ca}_{12}\text{Al}_{14}\text{O}_{32}\text{F}_2$ (C11A7f)	56
Gunnar Raade, Hans-Jørgen Berg	Powder X-ray diffraction data for innelite	62
Rashmi, D. K. Suri	X-ray powder diffraction study of CuInSeTe	65

INTERNATIONAL REPORTS


Regional Reports	69
Calendar of Meetings	72
Short Courses & Workshops	76
Cumulative Author Index	79

Enhance instrument **PERFORMANCE** with X-ray collimating lenses for parallel beam powder diffraction from **X-RAY OPTICAL SYSTEMS.**




x-ray collimator
(10x10 mm² output beam)
©XOS 1996

Parallel Beam Geometry



Bragg-Brentano Geometry



Outstanding Performance in **Stress and Texture Applications**


Elimination of all defocusing errors	No sample preparation needed
Large beam size	Improved particle statistics
Intensity gain up to 40 X	Improved counting statistics
2-D collimated quasi parallel beam	Increased orientation statistics
Constant peak profile and width over whole 2 Θ range	Ideally suited for industrial applications and full pattern analysis

Also available, X-ray focusing lenses for Micro X-Ray Fluorescence.
Call J.Phillip Bly, Sales and Application Engineer today.



30 Corporate Circle • Albany, NY 12203
Phone: 518.464.3334 • Fax: 518.464.3335
www.xos.com • email: info@xos.com

inel REAL TIME XRD



MPD Multi-Purpose X-ray Diffractometer

Versatile diffractometers designed to take advantage of the rapid, real time data collection our patented CPS X-ray detectors offer.


- *No scanning feature - acquire up to 120° 2 θ simultaneously
- Unique capillary devices for analysis of air sensitive materials
- Identify materials even if only micro amounts are available
- Parallel beam with germanium or multilayer mirror optics

for information on our complete product line please contact us

inel - Z.A - CD 405 - 45410 ARTENAY (FRANCE)
Tel. (33) 0 2 38 80 45 45 Fax. (33) 0 2 38 80 08 14
E.MAIL: inel@valcofim.fr-INTERNET:<http://www.valcofim.fr/inel>

inel Inc. P.O. Box 147, STRATHAM, NH 03885 (USA)
TEL. (603) 778-9161 FAX. (603) 778-9171
E-MAIL: inelinc@aol.com

The Most Important Materials In XRF Spectroscopy...Next To Your Samples



Binder & Grinding Aid

Roll Film

Plastic Vials & Balls

SPEX CertiPrep Disposables

Pellet Caps

XRF Sample Cups

Pre-Cleaned Bottles

Pre-Cut Film

Fusion Flux

For accurate and uniform analytical results, you need reliable and consistent sample preparation. SPEX CertiPrep has been supplying XRF spectroscopists with quality sample preparation and handling products for over forty years.

Whether you press powders in dies or fuse them into glass discs, or if you use sample cells to run liquids, powders, pastes, or other materials, we have the supplies and equipment to get the job done. We can help you choose the proper binder, grinding aid, or borate flux for your particular sample. SPEX CertiPrep also provides the most complete line of mills, presses, dies, and fusion fluxers.

Our commitment to Total Customer Satisfaction means: "Serving our customers promptly and courteously, with unsurpassed technical and

applications support, and top quality, reliable SPEX CertiPrep products." Call, FAX, or e-mail for your copy of our catalog, *The Handbook of Sample Preparation and Handling*. **1 800 LAB-SPEX**



203 Norcross Avenue ■ Metuchen, NJ 08840 USA ■ 732-549-7144
Fax 732-603-9647 ■ SamplePrep@spexcorp.com ■ <http://www.spexcorp.com>

MNI

Quantitative



Weight Fraction	Phase Name
37.39%	43-1484 Aluminum Oxide (Corundum, syn)
22.45%	05-0664 Zinc Oxide (Zincite, syn)
18.88%	47-1743 Calcium Carbonate (Calcite)
21.28%	21-1276 Titanium Oxide (Rutile, syn)

MNI:	0.42	Std. Deviation:	0.027	Density:	4.25
Volume Fraction:	0.20	Std. Deviation:	0.016	RIR:	3.40
Weight Fraction:	0.21	Std. Deviation:	0.017	Peaks used:	6

Figure of Merit: 0.040

Buttons: Peak details, Print, Close

Database Cards in Queue	
43-1484 Aluminum Oxide (Corundum, syn)	Save...
05-0664 Zinc Oxide (Zincite, syn)	
47-1743 Calcium Carbonate (Calcite)	Load...
21-1276 Titanium Oxide (Rutile, syn)	

Extract Database Card	Card: 21-1276	Extract
	Search String:	Search
0%		
<input checked="" type="radio"/> ICDD NBS AIDS83 Database	<input type="radio"/> Search-Match Database	Create PDF2 SubFile...
<input type="radio"/> Traces User Database		
Generate Patterns	Cancel	

Standardless Quantitative XRD Analysis using your ICDD PDF-2 Database

- Unlock new potential in your Database and extend its use from traditional Qualitative ID to new Quantitative Analysis
- The MNI method added to Tracesv5 allows you to identify and to analyse compounds in a mixture all in the one superbly comprehensive package
- No knowledge of crystallography or the "Rietveld" method needed
- Achieve accuracies as good as other pattern synthesis methods, automatically and in a fraction of the time
- Multi - line averaging compensates for preferred orientation
- Allow for non - diffracting or amorphous component in a mixture.

Please visit our Website for further details.
A Complete description with details and examples is available on our website.

- www.difftech.com.au -

Diffraction Technology Pty.Ltd

Contact: Rod Clapp

Phone: 61 2 62428233

Fax: 61 2 62428266

Email: difftech@difftech.com.au

P.O Box 444 Mitchell A.C.T 2911

Canberra, AUSTRALIA

United States: Jerry Krause
Analytical X-Ray Instrument Sales
Phone: 303 277 1569

United Kingdom: Arthur Bagley
Hiltonbrooks Ltd.

Phone: 01477534140


Fax: 01477535041

All Trademarks are the property of their respective owners


microSOURCE®

X - R A Y G E N E R A T O R S

it beats all the rest...



...so put it to the test



Contact us for examples of the excellent data we have collected and to learn how we can help you to **maximise your laboratory's potential.**

Demonstrations available on request

The Microsource® X-ray Generator gives you a complete system for less than \$75,000*

complete system

including tube mount, tube shield, X-ray optic (micromirror®) with built in alignment for X-ray focusing

compatibility

proven interface to different systems

performance

proven 10^{10} photons $s^{-1} mm^{-2}$ together with low divergence for a variety of applications

future proof

modular design, rapid servicing, upgradability

safety

built-in shutter and X-ray safety systems

economy

low power consumption, low cost of ownership, small footprint

For enquiries within USA, Canada & Mexico contact:
Dr D Keith Bowen, Bede Scientific Incorporated
14 Inverness Drive East, Suite G-104, Englewood, CO 80112, USA
Tel: 303.790.8647 • Fax: 303.790.8648 • Email: sales@bede.com
Website: www.bede.com

For the rest of the world contact:
Dr Graham Fraser, Bede Scientific Instruments Ltd.
Bowburn South Industrial Estate, Bowburn, Durham, DH6 5AD, UK
Tel: +44 (0) 191 377 2476 • Fax: +44 (0) 191 377 9952
Email: microsource@bede.co.uk • Website: www.bede.co.uk

*Price F.O.B. Durham, UK



microSOURCE® is a Division of Bede Scientific Instruments Ltd

Get specialty beryllium prototypes with production processing in mind.

Driven to perform? Come to us. We'll put the metal to the metal.

When your production schedule is short, you want to get it done right the first time. You can eliminate a major uncertainty when you entrust your custom beryllium fabrication, joining and coating to us. We'll help you shine.

Here's what we can do for you:

- Depend on us as the only fully integrated source for beryllium sheet and foil products. All critical operations are performed in-house for seamless excellence.
- We're your partner in joining and coating metals. When your project involves joining a metal to a metal, you can count on us.
- Come to us for all your UHV beryllium product needs: x-ray windows, chambers, beam pipes and more.
- Rely on our expert engineering. We like challenges. Taking your design concept into reality would be our pleasure.

Count on us for your tough jobs. Keep in mind that we're the only fully integrated source for beryllium products for the analytical, medical, and scientific industries. Challenge us with your "cutting edge" requirements. We look forward to hearing from you.

TEL: 510-623-1500 • FAX: 510-623-7600

E-Mail: Electrofusion@BrushWellman.com

44036 South Grimmer Boulevard • Fremont, California 94538 • USA



BRUSHWELLMAN

ELECTROFUSION PRODUCTS

We're the beryllium window folks.

2000 ICDD X-ray Clinics

Learn...
from a team of experts in the field

Need basic
theory?


Calibration
problems?

Learn tricks of
the trade

Learn...
techniques and procedures



X-RAY FLUORESCENCE SPECTROMETRY:
FUNDAMENTALS **1-5 May 2000**
ADVANCED METHODS **8-12 May 2000**

X-RAY POWDER DIFFRACTION:
FUNDAMENTALS **5-9 June 2000**
ADVANCED METHODS **12-16 June 2000**




International Centre for Diffraction Data
12 Campus Boulevard
Newtown Square, PA 19073-3273 U.S.A.

Tel: +610-325-9814 Fax: +610-325-9823
E-mail: clinics@icdd.com Web: www.icdd.com




PHASE ANALYSIS
65/45 MIX OF ALPHA & DELTA
POLYMORPHS OF MANNITOL


Kratos/Shimadzu XRD-6000



NON-AMBIENT STUDY
ASPIRIN OVER A TEMPERATURE
RANGE OF 40°C - 180°C STEPS



PERCENT CRYSTALLINITY
CALCULATED BY TOTAL
PATTERN FITTING



FULL PATTERN HALL METHOD
SHOWING VARIATION OF
CRYSTAL SIZE

X-Ray Diffraction is a unique analytical tool for solid state pharmaceutical products to unambiguously characterize their solid state nature. Many parameters required by the FDA, legal patent issues and drug performance are only accessible by utilizing X-Ray Diffraction.

X-Ray Diffraction is a direct result of crystal structures being present in the pharmaceutical under study. As such, the parameters typically associated with crystal structure can be simply assessed. For example, once an active drug has been isolated, an indexed X-Ray powder diffraction pattern is required to secure a patent and protect a company's investment. Additionally, analysis of the drug under most environmental conditions reveals the formation of any polymorphs which could adversely effect the drug's performance and toxicity.

With the introduction of new systems specifically designed for the pharmaceutical industry, Kratos and Shimadzu have made X-Ray Diffraction cost effective and easy to use.

**XRD-6000 Designed for
Pharmaceutical Research**

- ◆ Polymorph Screening
- ◆ Percentage Crystallinity
- ◆ Morphology
- ◆ Crystal Structure Analysis
- ◆ Excipient Quantification
- ◆ Non-Ambient & Humidity Variation

Call or Write Today for More Information:



Kratos Analytical, Inc.

Phone: (914)426-6700

E-Mail: info@kratos.com

Internet: <http://www.kratos.com>


**Call Today to Request Your Copy of
X-Ray Diffraction for
Solid State Pharmaceutical Products**

©2000 Kratos Analytical

SPEED IN POWDER DIFFRACTION

IMAGING PLATE GUINIER CAMERA 670*

featuring the Low Temperature Device 670.4: Temperature range: 300 K.. 8 K



* FAST, FLEXIBLE & FREE OF $K\alpha_2$



HUBER
X-RAY DIFFRACTION EQUIPMENT

HUBER Diffraktionstechnik GmbH
Sommerstrasse 4
D-83253 Rimsting
Germany

Tel.: +49 (0) 80 51 - 68 78 - 0
Fax: +49 (0) 80 51 - 68 78 - 10
e-Mail: info@xhuber.com
URL: www.xhuber.com



Max-Flux™ Powder Diffraction Upgrade



Benefits

- Reduction of sample displacement errors
- Enables testing of irregularly shaped samples
- Increased resolution
- Higher flux
- Lower background

Applications

- Thin Film Analysis
- High Resolution
- Phase Analysis
- Stress / Strain Measurements
- Texture Analysis
- Point-Defect Analysis
- And more...

 **Osmic** Inc.

An Advanced Materials and Technology Company

1788 Northwood Drive
Troy, Michigan 48084-5532, USA

tel: 248.362.1290

fax: 248.362.4043

www.osmic.com


Your Software For X-Ray Analysis

Now optimized for the new, enhanced PDF database with more than 42,000 additional calculated patterns



Software package for data acquisition and analysis:

Diffractometer control, plot modules, peak search and calibration, qualitative and quantitative phase analysis, profile fitting and pattern simulation, lattice parameter refinement, crystallite size and micro stress. Ready to plug to your PHILIPS and SIEMENS diffractometer hardware.



Stand alone package for qualitative and quantitative phase analysis.

Extremely fast search algorithms for ICDD-PDF database and user created databases. Very exact quantitative analysis using the WHOLE PATTERN method.



Extension package for most third party X-ray software. Graphic-analytical, profile-analytical and crystallographic options.


a. wassermann

röntgenanalytik.meßsysteme.software

P.O. BOX 2631, D-87416 Kempten, Germany, Tel. (49)0831-79515, FAX (49)0831-79930, E-mail: RMSKempten@aol.com

Are you working with current data?

Update to Release 1999 of the Powder Diffraction File™, the world's premier X-ray powder diffraction database.



- more than 120,000 patterns representing unique phases
- **new for '99** – 2,500 experimental patterns compiled through the ICDD's Grant-in-Aid program and literature searches
- **new for '99** – 4,220 patterns calculated from the ICSD database
- a four-tiered editorial process to assure the highest standard of accuracy and quality

International Centre for Diffraction Data

12 Campus Boulevard
Newtown Square, PA 19073-3273 U.S.A.

WEB SITE: www.icdd.com ♦ E-MAIL: info@icdd.com

SALES: 610-325-9810
TECHNICAL: 610-325-9814
FAX: 610-325-9823




The *Proceedings*
of the
Denver X-ray Conferences
FROM 1957-1995 ARE
NOW AVAILABLE ON
ONE CD-ROM!!

**ADVANCES IN X-RAY ANALYSIS
VOLUMES 1 THROUGH 39**

- ▼ Page images of all manuscripts found in Volumes 1-39
- ▼ Searchable cumulative table of contents compiled from Volumes 1-39
- ▼ Convenient navigation abilities
- ▼ Functions in Windows®, UNIX®, and Macintosh® environments

**ORDER YOUR COMPREHENSIVE
CD-ROM**

**ADVANCES IN X-RAY ANALYSIS
VOLUMES 1-39
FOR ONLY \$350**



International Centre for Diffraction Data
12 Campus Boulevard
Newtown Square, PA 19073-3273 U.S.A.
Sales: 610-325-9810 Fax: 610-325-9823

Internet: INFORMATION@ICDD.COM
Web Site: www.icdd.com

it's

SEARLE

the people, the complex challenges,
the stimulating environment. The opportunity to work
independently while contributing to a team.
The opportunity to make a difference - to leave your mark.
It's about food, health, hope and an exciting future
as we merge with Pharmacia and Upjohn.

**Research Investigator: Powder X-Ray
Diffraction & Elemental Analysis**

We are seeking a candidate to investigate, develop, evaluate and implement new methods, procedures and techniques to improve the scope and quality of Powder X-ray Diffraction and Elemental Analyses in support of R&D. This will include developing methods to quantitate polymorphs in organic compounds, establishing analytical procedures, and analyzing samples. You will also develop and maintain procedures for calibration, maintenance and operation, train and develop analysts on instrumentation, procedures and methods as well as interact with scientists and colleagues to solve problems and apply instrument techniques and methods. You will also explore other related analytical techniques for possible addition to the laboratory such as ICP, ICPMS, dedicated halogen analyzers, micro water determination etc.

Ph.D. in chemistry with an emphasis on PXD and 3-5+ years' experience with PXD of organic/pharmaceutical chemicals is essential. Knowledge of organic chemistry, crystallization, the study of polymorphs and a familiarity with elemental analysis and elemental analyzers for CHNS and halogens would be ideal but not necessary. **Ad Code: 00-0306**

Optical Spectroscopist

The selected candidate will investigate, develop, evaluate and implement new methods, procedures and techniques to improve the quality and scope of optical (IR, Raman, UV) spectroscopy and thermal analysis support to R&D, as well as identify instrument and/or computer interface problems, troubleshoot and affect solutions. Other responsibilities will be to develop, maintain and periodically review calibration and maintenance procedures, design experiments in IR, Raman, UV and thermal analyses, and train and develop analysts on instrumentation, procedures and methods.

Ph.D. in organic chemistry or equivalent, with 3-5+ years' experience in optical spectroscopy (Raman, IR, UV) and thermal analysis (DSC & TGA) is essential, combined with experience in interpreting Raman, IR, UV and TA data. Knowledge of Raman, IR, UV, thermal analysis instrumentation and laboratory practices are also required. Knowledge of basic electronics and computer technology and some supervisory experience desired. Skills in interpreting IR spectra ideal. **Ad Code: 99-2536**

In addition to an open environment, we offer a competitive salary and benefits package including tuition reimbursement and 401(k). Please submit your resume in scannable format to: monsanto@aoon-hros.com, Ad Code (selected from above) on your subject line and cover letter, or mail to: Monsanto (select Ad Code), PO Box 1262, Findlay, OH 45840. You will receive a postcard acknowledging receipt of your resume. EEO/AA Employer M/F/D/V. Please visit us on the web at www.searle-ican.com

I
can
Make A
Difference.

Publishing Agreement

1. It is our custom to request authors to vest the worldwide copyright of their papers with the JCPDS—International Centre for Diffraction Data for the full term of copyright and we would be grateful if you would confirm your acceptance of these terms by signing and returning the agreement below. We will not withhold permission for any reasonable request from you to publish any part of this paper in connection with any other work by you, provided the usual acknowledgements are given regarding copyright notice and reference to the original publication.
2. If it is appropriate, the author's employer may sign this agreement, and in any event the employer may reserve the right to use the paper internally or for promotional purposes only, by so indicating on this agreement. It is understood that proprietary rights other than copyright (including patent rights) are reserved.
3. If the author is a U. S. Government employee and this work was done in that capacity the assignment applies only to the extent allowable by U. S. law. If at least one co-author is *not* a government employee, said author should sign the agreement.
4. If the author is an employee of the British Government then HMSO will grant a non-exclusive license to publish this paper in the Journal in any form or media provided British crown copyright and user rights (including patent rights) are reserved.
5. The author warrants that the manuscript is the author's original work, and has not been published before. (If excerpts from copyrighted works are included, the author will obtain written permission from the copyright owners and show credit to the sources in the manuscript.) The author also warrants that the article contains no libelous or unlawful statements, and does not infringe on the rights of others.
6. If the work was prepared jointly the author agrees to inform co-authors of the terms of the agreement, and to sign on their behalf.

With the qualifications listed above I assign to JCPDS—International Centre for Diffraction Data, the copyright of my paper entitled

for publication in *Powder Diffraction*.

Signed

Date

If the paper is rejected, this assignment is null and void.