Part II: Total inward leakage tests. Int J Hyg Environ Health 2018;221: 977–984.

- 7. Chu DK, Akl EA, Duda S, *et al.* Physical distancing, face masks, and eye protection to prevent person-to-person transmission of SARS-CoV-2 and
- COVID-19: a systematic review and meta-analysis. *Lancet* 2020;395: 1973-1987.
- Advice on the use of masks in the context of COVID-19: interim guidance. World Health Organization website. https://apps.who.int/iris/handle/10665/ 332293. Published 2020. Accessed July 1, 2020.

Effect of vaporized hydrogen peroxide reprocessing on N95 respirators

Elena Beam MD¹ ⁽⁶⁾, Jeffrey C. Nesbitt CIH², Matthew D. Austin CIH² and Kannan Ramar MD³

¹Division of Infectious Diseases, Mayo Clinic, Rochester, Minnesota, ²Occupational Safety, Rochester, Minnesota and ³Division of Pulmonary and Critical Care Medicine, Center for Sleep Medicine, Mayo Clinic, Rochester, Minnesota

To the Editor—The high demand for personal protective equipment (PPE) during the coronavirus disease 2019 (COVID-19) pandemic has required reprocessing and reuse of N95 respirators to mitigate shortages. Data on the impacts of reprocessing and reuse on the physical integrity and continued performance of these PPE are sparse.

Our facility uses vaporized hydrogen peroxide (VHP) according to the Duke method as a strategy to reprocess respirators for reuse.¹ We conducted repeated quantitative fit testing of N95 respirators (model 3M 1870+, 3M, Maplewood, MN) by measuring the amount of leakage into the facepiece. We sought to better understand the impact of VHP reprocessing on reuse and extended reuse of respirators as well as the effect on the tight fit of the respirator in 2 experiments.

In our first experiment, 5 masks that were reprocessed with VHP were compared to 5 masks that were not treated with VHP. Quantitative fit testing was conducted using TSI PortaCount according to the Occupational Safety and Health Administration (OSHA) quantitative fit testing protocol. Safety staff conducted repeat use and fit tests in 10 cycles as follows: (1) Inspect the respirator. (2) Adjust and don the respirator. (3) Perform a quantitative respirator fit test. (4) Doff the respirator. And (5) flatten the respirator. Any observable failures and fit-test failures were documented after each cycle (Table 1). We found no fit-testing failure among the respirators treated with VHP, while the 5 control respirators had 4 total fit-test failures. Notably, 3 failures occurred after the ninth cycle of donning and doffing the respirators.

In our second experiment, to better understand the impact of VHP reprocessing on reuse and extended use (defined as repeat half-day use, or 4-hour shifts), we tested the same type of respirators (3M 1870) according to the following procedure: (1) The respirator was used for 4 hours. (2) The respirator was reprocessed with VHP. (3) The respirator was used for an additional 4 hours. (4) The respirator was reprocessed with VHP. And (5) a quantitative fit test was performed on the respirator. Respirators used in this experiment had already been used in hospital service by a single user for a clinical purpose. The duration of use was unknown.

All 5 respirators successfully passed the quantitative fit testing. The results from both experiments suggested that reprocessing with VHP allows for reuse and extended use of respirators. Our study has several limitations. These results might not be generalizable to other contexts. The quality and integrity of respirator used may differ by brand and style. Variability in the results may have resulted from the use of quantitative versus qualitative fit-testing techniques, which may be more prone to error due to the subjective nature of the test. Our evaluations were based on conclusions from quantitative testing.^{2.3} Finally, fit-testing expertise among the staff who volunteered to run the experiment varied; however, we expect our data to be reliable because the same staff performed multiple tests. Finally, we were not able to capture the duration of clinical use of the respirators in the hospital prior to the start of our second experiment, and we were unable to define the overall "use" in the project given this limitation.

Limited evidence is available regarding the use of VHP with extended reuse of N95 respirators.^{4,5} Similar to other studies, we found no detrimental effect of VHP processing on the ability of N95 respirators to pass fit testing. Our results suggest that VHP does not affect limited reuse and extended use of 3M-1870 respirators in the context of maintenance fit testing with repeated donning and doffing.

Cite this article: Beam E, et al. (2021). Effect of vaporized hydrogen peroxide reprocessing on N95 respirators. Infection Control & Hospital Epidemiology, 42: 907-908, https://doi.org/10.1017/ice.2020.371

© 2020 by The Society for Healthcare Epidemiology of America. All rights reserved. This is an Open Access article, distributed under the terms of the Creative Commons Attribution licence (http://creativecommons.org/licenses/by/4.0/), which permits unrestricted re-use, distribution, and reproduction in any medium, provided the original work is properly cited.

Author for correspondence: Elena Beam, E-mail: beam.elena@mayo.edu

Table 1. Respirator Fit Test Cycle Results

		Fit Test Cycle									
Respirator No.	1	2	3	4	5	6	7	8	9	10	
Control 1ª	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	
Control 2	Pass	Pass	Fail	Pass	Pass	Pass	Pass	Pass	Fail	Pass	
Control 3	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	
Control 4	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	
Control 5	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Fail	Fail	
VHP 1 ^b	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	
VHP 2	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	
VHP 3	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	
VHP 4	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	
VHP 5	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	Pass	

^aControl: 1870 respirator used without reprocessing.

^bVHP: 1870 respiratory used with vaporized hydrogen reprocessing.

Acknowledgments.

Financial support. No financial support was provided relevant to this article.

Conflicts of interest. All authors report no conflicts of interest relevant to this article.

References

- Schwartz A, Steigle M, Greeson, N, et al. Decontamination and reuse of N95 respirators with hydrogen peroxide vapor 2 to address worldwide personal protective equipment shortages during the 3 SARS-CoV-2 (COVID-19) pandemic. Appl Biosafety 2020 [Epub ahead of print]. doi: 10.1177/ 1535676020919932.
- 2. Hon CY, Danyluk Q, Bryce E, *et al.* Comparison of qualitative and quantitative fit-testing results for three commonly used respirators in the healthcare sector. *J Occup Environ Hyg* 2017;14:175–179.
- Occupational safety and health standards: occupational health and environmental control (standard no. 1910.134). Occupational Safety and Health Administration website. https://www.osha.gov/laws-regs/regulations/standardnumber/1910/1910.134AppA. Published 1970. Accessed July 23, 2020.
- Fischer RJ, Morris DH, van Doremalen N, et al. Assessment of N95 respirator decontamination and re-use for SARS-CoV-2. medRxiv 2020 Apr 24. doi: 10.1101/2020.04.11.20062018.
- Degesys NF, Wang, RC, Kwon, E, *et al.* Correlation between N95 extended use and reuse and fit failure in an emergency department. *JAMA* 2020;324: 94–96.

RE: Universal SARS-CoV-2 testing on admission to the labor and delivery unit: Low prevalence among asymptomatic obstetric patients

Sean Cronin MD, Megan Piacquadio DO, Katelyn Brendel DO, Aden Goldberg, Marco Goldberg, Chase White MD, David Jaspan DO, and Jay Goldberg MD 💿

Department of Obstetrics and Gynecology, Einstein Medical Center Philadelphia, Philadelphia, Pennsylvania

To the Editor—In their recent publication, Goldfarb et al¹ reported a low prevalence of coronavirus disease 2019 (COVID-19), 1.5%, among asymptomatic pregnant women in Boston presenting for admission to labor and delivery between April 18, 2020, and

Author for correspondence: Jay Goldberg, E-mail: jaygoldbergmd@yahoo.com

Cite this article: Cronin S, et al. (2021). RE: Universal SARS-CoV-2 testing on admission to the labor and delivery unit: Low prevalence among asymptomatic obstetric patients. Infection Control & Hospital Epidemiology, 42: 908–909, https://doi.org/10.1017/ice.2020.382

May 5, 2020. Noting that their rate was substantially lower than that reported in New York City, the authors theorized that it might be due to their patients (1) being tested >30 days after physical distancing orders were in place; (2) the population density of Boston being less than New York City; and (3) New York women underreporting symptoms due to New York hospitals banning support people from labor and delivery.¹

Studying similar universal screening in pregnant women presenting to labor and delivery at Einstein Medical Center Philadelphia during the same time frame as the Boston study,

© 2020 by The Society for Healthcare Epidemiology of America. All rights reserved. This is an Open Access article, distributed under the terms of the Creative Commons Attribution licence (http://creativecommons.org/licenses/by/4.0/), which permits unrestricted re-use, distribution, and reproduction in any medium, provided the original work is properly cited.