Advertisers in This Is	SSUC Page No.
* Aldrich Chemical Company, Inc.	54
* American Institute of Physics	64
* Amptek Inc.	33
Annual Reviews	73
Arizona State University	7
Axic, Inc.	58
* Bio-Logic Science Instruments	20
* Cambridge University Press	101
* Cerac Inc.	76
Chemat Technology, Inc.	47
* Digital Instruments	10
* Edax International Inc.	14, 15
* Elsevier Science	84
* ESM Software	53
* Gordon & Breach	82
* Hinds Instruments, Inc.	7
* High Voltage Engineering	inside front cover
* Huntington Laboratories	outside back cover
* JEOL USA, Inc.	16
* Lambda Technologies Inc.	26
* MDC Vacuum Products Corporation	25
The Mellen Company, Inc.	73
* MMR Technologies, Inc.	20
* Nanopowder Enterprises Inc.	11
* National Electrostatics Corporation	9
New Focus, Inc.	inside back cover
* NFT—Nanofilm Technology	64
* Omicron Associates	31
* Osmic, Inc.	24
Outokumpu Research Oy	13
Oxford Cryosystems	17
* Oxford University Press	102
* Park Scientific Instruments	1
* Philips Analytical X-Ray	23
* Plenum Publishing Corporation	74
* Praxair Specialty Ceramics	21
* Princeton Gamma-Tech Inc.	18
* Quad Group	27
* Quantum Design	8
* Solartron	2
* Specs USA, Inc.	78
* Springer-Verlag	77, 100
* Structured Materials Industries, Inc.	11
* Sycon Instruments	6
* VAT, Inc.	12
* Virginia Semiconductor, Inc.	32, 79
VLSI Standards Inc.	78
* Voltaix, Inc.	28


* Please visit us at the Exhibit, December 2–4, during the 1997 MRS Fall Meeting in Boston.

For free information about the products and services offered in this issue, fill out and mail the Reader Service Card, or FAX it to 312-922-3165.

New Solutions


Materials and Surface Analysis

SNOM - The New Dimension in Optical Resolution



ESCAPROBE - Surface Analysis and Materials Research System

• Uncompromised ESCA Performance



Headquarters: OMICRON VAKUUMPHYSIK GmbH

Idsteiner Str. 78 D-65232 Taunusstein, Germany Tel. +49 (0) 6128 987-0 Fax +49 (0) 6128 987-185

USA: OMICRON Associates 1226 Stoltz Road Bethel Park, PA 15102

Tel. (412) 831-2262


Fax (412) 831-9828

INSTRUMENTS FOR SURFACE

http://www.omicron-instruments.com/

"But still try—for who knows what is possible?


- Michael Faraday



10-MICRON THICK SILICON WAFER



4000-MICRON THICK SILICON WAFER



MICROMACHINED SILICON WAFER



The ubiquitous nature of single crystal silicon provides for application possibilities which go far beyond those defined by the "traditional" microelectronics industry. Virginia Semiconductor, Inc. considers silicon to be an **ENGINEERED MATERIAL** ideal for a host of applications that call for

- fatigue and chemical resistance
- machinability
- mechanical robustness
- · thermal stability
- · and electrical conductivity.

When it comes to preparing engineered silicon products, we have frequently stated (with tongue in check), "if we can't make it, you don't need it!"


In the final analysis, we are most eager to **IMAGINE**, to be **CHALLENGED**, and to TRY. At Virginia Semiconductor, Inc., we think of the possibilities—not the limitations. Those who know us now expect nothing less; why shouldn't it be that way in our service to you?



Virginia Semiconductor Inc.

1501 Powhatan Street Fredericksburg, VA 22401

Phone: (540) 373-2900 Fax: (540) 371-0371



XR-100CR

X-RAY & GAMMA RAY




APPLICATIONS

- X-Ray Detection
- FLUORESCENCE
- ENVIRONMENTAL MONITORS
- NUCLEAR MEDICINE
- PORTABLE INSTRUMENTS
- SPACE & ASTRONOMY
- TEACHING AND R&D

FEATURES

- . SOLID STATE DESIGN
- THERMOELECTRIC COOLER
- HERMETIC PACKAGE
- WIDE DETECTION RANGE
- BERYLLIUM WINDOW
- Easy to Operate



FIRST X-RAY SPECTRUM FROM MARS

4 Energy (keV)

For its unique design and reliability, this detector was selected for the Mars Pathfinder Mission and placed on the Rover arm to perform rock and soil analysis via x-ray fluorescence techniques.


10

NEW!!! 196 eV Resolution **XR-100CR** With Si-PIN for X-Ray Detection


⁵⁵Fe Spectrum (20 μs shaping) C 0 U 196 eV FWHM Ν Т S 6.4 keV ENERGY (keV)

XR-100T-CZT

With Cadmium Zinc Telluride (CZT) for γ-Ray Detection



Models XR-100CR and XR-100T-CZT are high performance X-Ray and γ-Ray detectors mounted on a thermoelectric cooler together with the input FET to the Preamplifier. Monitored by an integrated circuit, these components are kept at -30°C and are enclosed in a hermetic TO-8 package with a vacuum tight, light tight Beryllium window. Power and signal processing to the detector is provided by the PX2T in order to ensure quick, stable operation in less than one minute from power turn-on. The output pulse produced by the PX2T can be connected directly to the input of a Multichannel Analyzer (MCA). For optimum portability and versatility, use the Amptek MCA8000A "Pocket MCA."



AMPTEK INC. 6 DE ANGELO DRIVE, BEDFORD, MA 01730-2204 U.S.A. Tel: +1 (781) 275-2242 Fax: +1 (781) 275-3470 e-mail: sales@amptek.com

Visit Us Now www.amptek.com

