

The Marconi Pure Tone Audiometer is now available. It has been developed for the Ministry of Supply to meet the specification proposed by the Committee on Electro-Acoustics appointed by the Medical Research Council. Each instrument is standardized against N.P.L. determinations so that, from a given patient, substantially comparable results can be obtained with any Marconi Audiometer. The patient, who can be orally instructed through a crystal microphone and high fidelity ear-phones, communicates with the operator via a built-in signal lamp. The sound output over the normal hearing range is automatically corrected for characteristic losses.

THE PURE TONE

Audiometer

Type TF 895

For full particulars of this efficient aid to audiometry please write:

MARCONI INSTRUMENTS LTD.

ST. ALBANS, HERTFORDSHIRE . Telephone: St. Albans 6161/5

Northern Office: 30 Albion Street, Hull. Western Office: 76 Portview Road, Avonmouth. Midland Office: 19 The Parade, Leamington Spa. Southern Office and Showrooms: 109 Eaton Square, London, S.W.1.



