# The Journal of Laryngology and Otology

G. H. BATEMAN

ASSISTANT EDITOR LIONEL TAYLOR

#### Contents


The scope and application of current audiometric		
tests	Gunnar Lidén	507
Auditory analysis in some common hearing prob-		
lems	T. J. WILMOT	521
Histopathological study of the middle-ear cleft .	L. H. HIRANANDANI	
	and C. K. DESHPANDE	529
Congenital atresia of the ear-a review of the		
surgical findings in 83 cases	N. W. GILL	55I
Primary myringotomy for secretory otitis media in		
children	R. J. BENNETT	
	and A. N. CHARRABORTY	589
Dissection tonsillectomy-complications and		
follow-up	KURT BRESON	
	and JAAP DIEPEVEEN	DOI
Clinical records—		
The external carotid artery (a perfusion technique		
used in 30 patients with head and neck cancer)	WILLIAM S. LUND	609
Infectious mononucleosis, deafness and facial		
paralysis	LIONEL TAYLOR	
	and GERALD PARSONS-SMITH	613
Respiratory obstruction in glandular fever	M. D. LEE	617
Chondroma of larynx	K. D. GUPTA	623
Fibrosarcoma of the ethmoid	U. PRASAD	
	and J. K. KANJILAL	627
Societies' proceedings		633

### Headley Brothers Ltd

Ashford, Kent

General notes

# **Announcing...** an Audiometer of exceptional distinction



# THE AMPLIVOX MODEL 103 DIAGNOSTIC AUDIOMETER

to ISO Specification R389:1964


A sophisticated design for the comprehensive evaluation of hearing acuity and the differential diagnosis of hearing impairments

- \* Easy-to-use controls for all tests from simple clinical to advanced research levels
- \* Narrow band masking for pure tone, wide band white noise for speech
- \* Unique masking threshold indicator
- \* Plug-in Adaptor for SISI and Lüscher Tests
- \* Unique Patient's Signal assists subject's concentration
- \* Photo-electric switching for noise-free operation
- \* Twin channel facilities with multiple switching for dozens of tests
- \* Plug-in circuit boards for ease of service
- \* Output for free field audiometry

Please write for technical literature

# **AMPLIVOX**

Hearing Centres, Ltd Beresford Avenue, Wembley, Middx. Tel: 01 - 902 8991



### WHEN A HEARING AID IS RECOMMENDED

# REFER YOUR PATIENTS TO THE EXPERTS—

### **AMPLIVOX**

#### **Hearing Aids**

A wide selection of hearing aids always available from the world's manufacturers. We are the exclusive agent in the United Kingdom and Eire for Dahlberg, U.S.A., in addition to Amplivox.

#### **Ear Defenders**

All types of ear defenders available at all branches.

#### **Medical Acoustic Division**

Audiometers, Acoustic Booths, Group and Speech Training Hearing Aids. We are the exclusive agent for the United Kingdom and Eire for Tracor, U.S.A.

FOR YOUR PATIENTS' GREATEST SATISFACTION SPECIFY AMPLIVOX HEARING CENTRES, THE ONLY COMPANY IN THE WORLD MANUFACTURING AND SUPPLYING A COMPLETE RANGE OF EQUIPMENT FOR THE PREVENTION, TESTING AND ALLEVIATION OF DEAFNESS

#### AMPLIVOX HEARING CENTRES 80 NEW BOND STREET, LONDON, W.1. TEL.: 01-493 9888

#### Branches at:

London · Birmingham · Bournemouth · Bristol · Cardiff · Dublin Edinburgh · Glasgow · Leeds · Liverpool · Manchester · Newcastle · Wembley

Specifically designed for you and your patients. Precision-engineered, functionally-built equipment for the Ear, Nose and Throat and other Specialists.

S.M.R. keeps you and your patients' comfort and satisfaction in mind. Whatever your needs . . . whether they be for one unit or a completely equipped treatment room, you are assured of the finest . . . a room that reflects professional good taste.



Picture shows what can be done with small treatment room areas. The rooms are complete from the compact treatment unit to the physician's stool. Yet, this equipment fits easily and functionally into a small area. This is just one of many room arrangements available from S.M.R.




#### STAINLESS Steel Polish

SMR Stainless
Steel Polish
cleans and preserves with
superior results
—ends fingerprints,
watermarks and
surface stains.
22 oz. economy

size . . . . \$2.50

#### SURGICAL MECHANICAL RESEARCH INC.

P.O. Box 1185, Newport Beach, Calif. 92663 FACTORY DIRECT SALES



Please send for our new catalog Direct factory sales, F.O.B. Newport Beach, Calif.

# Refer your patients to Inglieur patients to Inglieur patients to

# Ingrams have the different makes to choose from

LONDON 2 Shepherd Street, London, W.1. Hyde Park 9041

(2 minutes from Green Park)

BIRMINGHAM 20 Auchinleck Sq, Broad St, Birmingham, 15. Midland 4314

(In the Five Ways Shopping Centre)

BRIGHTON 5 Duke Street, Brighton, 1, Brighton 28699

(First turning off West Street)

CROYDON 11 Essex House, Co'lege Road, Croydon. Croydon 4663

(Opposite the Technical College)

MANCHESTER Mansfield Chambers (Rooms 21 & 22), Blackfriars 4104

17 St. Anne's Square, Manchester, 2. (STD Code (Off Deansgate) 061-834 4104)

#### ...and in SCOTLAND it's BONOCHORD-INGRAMS...

ABERDEEN 461 Union Street, Aberdeen. Aberdeen 29632

(Near Holborn Junction)

EDINBURGH 17 Stafford Street, Edinburgh. Caledonian 1942


(Corner of Alva Street)

GLASGOW 38 Bath Street, Glasgow, C.2. Douglas 7710

(Next door to Glasgow Corporation Transport Office)

Domiciliary visits in EVERY part of the country.

INGRAMS FOR THE LARGER RANGE OF HEARING AIDS—
Subscribers to the Code of Ethics of the Hearing Aid Industry Association.



#### PETERS CLINIC AUDIOMETERS AND BÉKÉSY


ARE THE IDEAL TOOLS FOR MORE ADVANCED AUDIOMETRY

#### The Audiometers offer:

- Two continuously variable attenuators accurate to I db.
- SISI Test (built-in) and automatic alternate loudness balance switching.
- Maximum attenuation of 120 db for loudness discomfort level tests.
- Narrow band SYNCHRONOUS masking (a new semi-automatic masking technique).
- Full speech test facilities with white noise masking.
- Plug-in circuit boards and operator output checks.
- Conversion to a true BÉKÉSY AUDIOMETER at will.
- Extra features not listed are offered by the AP/6 compared to the AP/5 AUDIOMETER.

#### write for full details to:

# ALFRED PETERS AND SONS LTD GELL STREET SHEFFIELD 3



Clinicians have frequently testified to the success of Orbenin in the treatment of Gram-positive infections. World Literature now contains over 400 references to Orbenin.

Pneumonia 🗌 Lung Abscess 🔲 Tonsillitis 🔝 Pharyngitis 🗀 Acute bronchitis	
Septicaemia	-
Boils 🗌 Carbuncles 🔲 Infected dermatoses	-
Osteomyelitis 🗌 Osteitis	
Post-operative wound infections 🗌 Burns 🔲 Skin graft protection	
Acute endocarditis	
Staphylococcal enterocolitis  Staphylococcal urinary tract infections Staphylococcal meningitis	

#### Dosage

Adults: 500mg (2 capsules) four times daily, by mouth. 250mg (1 vial) by intramuscular injection, four times daily.

Children: Under 2 years—quarter adult dose. 2-10 years—half adult dose.

In severe infections the dosage may be increased.

#### Contra-indications

Orbenin should not be given to patients with a penicillin allergy or administered by subconjunctival injection.

#### Side-effects

As with other penicillins.

Additional information is available on request.



Orbenin (Regd.) (cloxacillin sodium BP) is a product of British research at **Beecham Research Laboratories** Brentford, England. originators of the new penicillins.



#### THE LARYNGOSCOPE

A Monthly Journal devoted to the disease of EAR. NOSE AND THROAT

Official organ for the American Laryngological Rhinological and Otological Society

Price \$18.00 per year Canada \$19.00 per year Foreign \$19.00 per year

ESTABLISHED 1896

THEODORE E. WALSH, M.D. EDITOR

517 SOUTH EUCLID AVENUE SAINT LOUIS MO. 63110.



#### P. C. WERTH LTD. Contractors to H.M. Government and Local Authorities

Specialist Supply, Maintenance, Calibration Services for Audiological Test and Teaching Equipment

STRATFORD PLACE, OXFORD ST., LONDON W1N ODH

Telephone: 01-629 5559

Clinical, Diagnostic, Screening and Bekesy Audiometers. Acoustic Impedance Apparatus. Stereo Auditory Trainers. Visible Speech Apparatus. Artificial Mastoid, etc.

Our latest model with many new features:

New audiogram light indi-


Masking to 100dB intensity for pure tone and speech. Fowler, Sisi and Luscher


Automatic pulsing with 4 time intervals.

Complete speech facilities with 4 built-in weighting curves.

Simplicity of operation.

Probably the most keenly priced twin-channel audiometer currently available.





# THE DROP

# THAT SOOTHES AND CLEARS THE INFECTED EAR

'Otosporin' rapidly destroys virtually *all* bacteria commonly found in ear infections, reduces irritation, inflammation and pain. 'Otosporin' brand Drops contain polymyxin B sulphate, neomycin sulphate and hydrocortisone in a bland suspension designed to facilitate penetration. Issued in bottles of 5 ml.

### **'OTOSPORIN'**

Full information is available on request.



Burroughs Wellcome & Co. London
The Wellcome Foundation Ltd.)

#### eliminate the cross infection risk



# Opulets sterile, single dose eye and ear preparations

To overcome the risk of contamination in eye and ear infections, Opulets have been developed as sterile, single dose eye and ear drop tubes and ointment capsules. Each is packed in a film sachet and sterilized by gamma radiation, ensuring effective sterilization with no loss of activity—the strict criteria of efficacy and safety being met in full.

After each application, these single dose units are discarded, eliminating cross infection risks – each treatment administered by either medical staff or the patient is thus under the most sterile conditions possible in practical medicine.

Opulets provide the lowest unit price available compared with all other preparations of a similar nature. Presentations available at present are: Atropine Eye Ointment Atropine Sulphate BP 0.5% w/w. Chloramphenicol Eye and Ear Drops Chloramphenicol BP 0.5% w/v. Chloramphenicol Eye and Ear Ointment Chloramphenicol BP 1.0% w/w. Fluorescein Eye Drops Fluorescein Sodium BP 1% w/v Hydrocortisone/Neomycin Eye and Ear Ointment Hydrocortisone Acetate BP 1.5% w/w Noemycin BP 0.5% w/w. Physostigmine Eye Drops Physostigmine (Eserine) Sulphate BPC 0.5% w/v. Sulphacetamide Eye Drops Sulphacetamide Sodium BP 10% w/v.

Sulphacetamide Eye Ointment

Sulphacetamide Sodium BP 6% w/w.

full descriptive technical literature is available on request Pharmax Limited Thames Road Crayford Kent

